

# Performance of PLASTEM<sup>®</sup> as a Consistent Culture Supplement for Cell Therapy Manufacture

Irene Oliver-Vila<sup>1</sup>, Francisco J. Belda<sup>2</sup>, Eduardo Sesma<sup>3</sup>, Anna Seriola<sup>1</sup>, Samuel Ojosnegros<sup>1</sup>

<sup>1</sup> IVFTECH Parc Científic de Barcelona, Barcelona, Spain; <sup>2</sup> BioSupplies Research & Development Department, Grifols, Sant Cugat del Vallès, Barcelona, Spain; <sup>3</sup> IBEC, Institute for Bioengineering of Catalonia, Barcelona, Spain

## INTRODUCTION

- Culture media supplements such as human platelet lysate (hPL), fetal bovine serum (FBS), and male AB serum (hSerAB) are essential raw materials for cell therapy manufacture.
- Supplements must be of high quality and compliant with regulatory standards, but they also present challenges due to safety concerns, potential complexity, individual batch validation and limited batch sizes<sup>1</sup>.
- PLASTEM<sup>®</sup> (Grifols) is a consistent culture media supplement manufactured by fractionation of human plasma from >1000 donors, with a robust safety profile, consistency and scalable production. PLASTEM<sup>®</sup> combined with fibroblast growth factor (FGF) is an effective media supplement for human mesenchymal stem cell culture<sup>2</sup>.

## AIM

To evaluate the performance of PLASTEM<sup>®</sup> as a culture media supplement for cell therapy manufacture.

## METHODS

Human mesenchymal stromal cells (hMSCs) were cultured in DMEM supplemented with either FBS 10%, PLASTEM<sup>®</sup> 10% + hPL 0.5%, or PLASTEM<sup>®</sup> 10% + FGF.

Human Peripheral Blood Pan-T Cells (T cells) were cultured with RPMI 1640 supplemented with either, FBS 10%, PLASTEM<sup>®</sup> 15%, or hSerAB 5%.

**Cell density** was determined by Trypan blue exclusion assay for hMSCs and T cells. Expansion factor was calculated as final/initial seeding concentration.


**Multipotentiality assays.** Osteogenesis (Alizarin Red S staining), chondrogenesis (dark-blue copper-containing dye Alcian Blue) and adipogenesis (Oil Red staining) were determined in hMSCs supplemented with FBS or PLASTEM<sup>®</sup> + hPL 0.5% after 21 days.

**Immunomodulation assays.** hMSCs supplemented with FBS, or with PLASTEM<sup>®</sup> 10% + hPL 0.5% were expanded for 21 days. Peripheral blood mononuclear cells (PBMCs) were stained with carboxyfluorescein succinimidyl ester, co-cultured with hMSCs and activated with phytohemagglutinin (PHA). The proliferation percentage was determined by flow cytometry.

**Immunophenotyping.** Expression of cell surface markers was determined using flow cytometry in cells supplemented with PLASTEM<sup>®</sup> 10% + hPL 0.5%. hMSCs were stained with anti-CD90, anti-CD73, anti-CD105, anti-CD45, anti-CD14, anti-CD19, anti-CD34, and anti-HLA-DR. T cells were stained with anti-CD3, anti-CD4, anti-CCR7, anti-CD45, anti-CD45RA and anti-CD8.


## RESULTS

- In hMSCs, PLASTEM<sup>®</sup> 10% + FGF resulted in an efficient **cell proliferation** (71±9.5%) after 13-day compared with FBS (52.4±1.36%) or PLASTEM<sup>®</sup>+ hPL 0.5% (39.2±4.6%) (Figure 1A). T cells grew at similar rate in both media (fold expansion of cells after 12-day culture; FBS (54.2±7.5%), hSerAB 5% (49.2±18.6%), and PLASTEM<sup>®</sup> 15% (54.9±4.7%) (Figure 1B).




**FIGURE 1.** Determination of cell proliferation in (A) umbilical cord human mesenchymal stromal cells (UC-hMSCs) supplemented with FBS 10%, PLASTEM<sup>®</sup> 10% + hPL 0.5% or PLASTEM<sup>®</sup> 10% + FGF and (B) Human Peripheral Blood Pan-T Cells (T cells) supplemented with FBS, PLASTEM<sup>®</sup> 15% or hSerAB 5%. Results expressed as fold expansion between final and initial cell density (n=3).

- hMSCs cultured in PLASTEM<sup>®</sup> + hPL 0.5% supplemented media retained its phenotypic characteristics: adherence to plastic, the **multipotentiality** into osteoblasts, adipocytes or chondrocytes (Figure 2A) and the expression of **specific immune markers** CD90, CD73 and CD105 (99.9-100% expression) (Figure 2B). Negative markers for hMSCs were CD45, CD14, CD19, CD34 and HLA-DR (<5% expression).



**FIGURE 2.** (A) Differentiation of umbilical cord human mesenchymal stromal cells (UC-hMSCs) into the osteogenic, chondrogenic and adipogenic lineages. (B) Immunophenotype flow cytometry histograms for UC-hMSCs supplemented with PLASTEM<sup>®</sup> 10% + hPL 0.5% determining the expression of CD markers CD90, CD73, and CD105.

- In **immunomodulation assays**, hMSCs supplemented with PLASTEM<sup>®</sup> 10%+ hPL 0.5% reduced PBMCs absolute proliferation percentage (31.5%, Promocell; 37.6%, Creative Bioarray) compared with PBMCs alone (82.2%) (Figure 3).



**FIGURE 3.** Immunomodulation assay by co-culture of PBMCs and UC-hMSCs. Absolute proliferation of PHA-activated PBMC alone or in a co-culture with UC-hMSCs from two suppliers (Promocell and Creative Bioarray) and supplemented with FBS or PLASTEM<sup>®</sup> 10% + hPL 0.5%.

- The **immunophenotype** of T cells supplemented with PLASTEM<sup>®</sup> 15% showed an increased proportion of CD4+ T helper cells compared with FBS or hSerAB (Table 1).

**Table 1.** Immunophenotyping of T cells cultured with different media supplements are determined by cell surface expression of CD markers

Cell population	CD Surface marker	% of T cells, range, n=2		
		FBS	hSerAB	PLASTEM 15%
T- helper cells (CD4+)	CD45 <sup>+</sup> CD3 <sup>+</sup> CD4 <sup>+</sup>	60.1 - 60.3	46.1 - 46.8	65.4 - 80.0
Cytotoxic T cells (CD8+)	CD45 <sup>+</sup> CD3 <sup>+</sup> CD8 <sup>+</sup>	9.6 - 37.3	25.4 - 52.1	3.47 - 16.1
Naïve T cells	CD45 <sup>+</sup> CD3 <sup>+</sup> CD4 <sup>+</sup> CCR7 <sup>+</sup> CD45RA <sup>+</sup>	3.4 - 3.7	1.1 - 3.1	1.1 - 3.1
	CD45 <sup>+</sup> CD3 <sup>+</sup> CD8 <sup>+</sup> CCR7 <sup>+</sup> CD45RA <sup>+</sup>	7.9 - 17.9	1.0 - 4.1	2.3 - 12
Effector T cells	CD45 <sup>+</sup> CD3 <sup>+</sup> CD4 <sup>+</sup> CCR7 <sup>+</sup> CD45RA <sup>+</sup>	57.3 - 60	60.1 - 68.9	73.9 - 76.6
	CD45 <sup>+</sup> CD3 <sup>+</sup> CD8 <sup>+</sup> CCR7 <sup>+</sup> CD45RA <sup>+</sup>	77.9 - 88.5	92.2 - 95.8	87.8 - 94.2
Effector memory T cells	CD45 <sup>+</sup> CD3 <sup>+</sup> CD4 <sup>+</sup> CCR7 <sup>+</sup> CD45RA <sup>-</sup>	34.4 - 35.9	27 - 38.2	21.0 - 21.7
	CD45 <sup>+</sup> CD3 <sup>+</sup> CD8 <sup>+</sup> CCR7 <sup>+</sup> CD45RA <sup>-</sup>	2.5 - 3.8	3.1 - 3.5	0.2 - 3.3
Central memory T cells	CD45 <sup>+</sup> CD3 <sup>+</sup> CD4 <sup>+</sup> CCR7 <sup>+</sup> CD45RA <sup>+</sup>	1.8 - 3.6	0.7 - 1.0	1.2 - 1.3
	CD45 <sup>+</sup> CD3 <sup>+</sup> CD8 <sup>+</sup> CCR7 <sup>+</sup> CD45RA <sup>+</sup>	0.3 - 1.1	0.2 - 0.2	0 - 0.4

## REFERENCES

- Trento, C. et al. *Biol. Blood Marrow Transplant.* 2018;24:2365–70
- Blazquez-Prunera, A. et al. *Stem Cell Res Ther.* 2017;8(1):103

## DISCLOSURES

Francisco J. Belda is a full-time employee of Grifols, a manufacturer of PLASTEM<sup>®</sup>.

## CONCLUSIONS

PLASTEM<sup>®</sup> demonstrated suitability and consistency for cell therapy manufacturing process and it was proved as a competitive and robust supplement for cell culture.

The culture media supplement will help to streamline regulatory submissions, enabling the reproducible generation of potent quality-assured cell therapy.